Manual de usuario • User guidelines • Manual d'usuari
Manuel d'utilisation • Gebrauchsanweisung • Manuale d'uso
Gebruiksaanwijzing • Instrukcja obsługi • Manual de utilização
Brugsanvisning • Användarmanual • Brukerhåndbok
Käyttöopas • Notendahandbók • Lietotošanas pamācība
Ghidul utilizatorului • Priročnik za uporabo • Εγχειρίδιο λειτουργίας

HANDBIKE BATEC ELECTRIC

- AIR
- MINI 2 PLUS
- ELECTRIC 2 PLUS (RAPID 2 PLUS)
- SCRAMBLER 2 PLUS
- RAPID SPORT
- SCRAMBLER SPORT

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HANDBIKE BATEC ELECTRIC

- AIR
- MINI 2 PLUS
- ELECTRIC 2 PLUS (RAPID 2 PLUS)
- SCRAMBLER 2 PLUS
- RAPID SPORT
- SCRAMBLER SPORT

ÍNDICE

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CONTENTS

1. Product description, intended purpose, indications and contraindications	5
2. Safety measures and general precautions	7
2.1. Before using your HANDBIKE BATEC ELECTRIC for the first time	7
2.2. During use	8
2.3. General precautions	11
2.4. Lifespan	12
3. Models: characteristics and parts	13
3.1. HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT	13
4. Operating the BATEC HANDBIKE / user guide	18
4.1. Delivery	18
4.2. Unpacking	18
4.3. Elevation, transport and storage	18
4.4. Attaching the HANDBIKE BATEC ELECTRIC	19
4.5. Detaching the HANDBIKE BATEC ELECTRIC	19
4.6. Folding and unfolding the HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS	20
4.7. Using the HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT controls	21
4.8. Using the HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS/ RAPID SPORT / SCRAMBLER SPORT QUAD controls	24
4.9. Operating the HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS HEMI controls	27
4.10. Anti-Slip System	30
4.11. Using the Cruise Control System	31
4.12. Battery	32
4.13. Weights	37
5. Installation and adjustments	38
5.1. Installation of the HANDBIKE BATEC ELECTRIC on a wheelchair	38
5.2. User adjustments	38
6. Maintenance and repairs	43
7. Label	46

8. Warranty	47
8.1. Warranty period	47
8.2. Warranty conditions	47
8.3. Warranty card	48
9. Revisions list	49

INTENDED PURPOSE, PRODUCT DESCRIPTION, INTENDED USE , INDICATIONS, CONTRAINDICATIONS AND SPEED

INTENDED PURPOSE

HANDBIKE BATEC ELECTRIC products are accessories for medical devices that are intended to facilitate the personal mobility of users with total or partial disability in the lower limbs. They are designed to be attached to a manual wheelchair, allowing them to be moved with the assistance of an electric motor located on the front wheel.

PRODUCT DESCRIPTION

HANDBIKE BATEC ELECTRIC has CE marking as a class I medical device according to the requirements of Regulation (EU) 2017/745 concerning medical devices.

All the HANDBIKE BATEC ELECTRIC models operate in a similar manner, the main differences being the wheel size, top speed and features. See Section 3 (MODELS).

First, the handbike must be attached to a wheelchair. Once the two parts have been joined, we have a three-wheeled vehicle with the seat on the wheelchair and transmission and traction on the front wheel. To operate the handbike, just activate its transmission by pedaling with your hands or pressing on the accelerator. This will provide traction to the front wheel.

INTENDED USE

This product aims to provide wheelchair users (people with a disability that makes it impossible for them to walk) with practical and functional products that are useful both in their daily life and for doing recreational activities and sport.

Thanks to the product configuration, users can do longer and more complicated trips than with a wheelchair alone, boosting independence and autonomy. The Sport model enables higher speeds when used on closed tracks or other controlled environments.

INDICATIONS

The HANDBIKE BATEC ELECTRIC is indicated for users with the following characteristics:

- With a disability affecting the lower limbs, who use a wheelchair and whose upper limbs are not affected and have full functionality of their hands. For these users, models with either the standard or PARA set-up are indicated.
- With a disability affecting the lower limbs, who use a wheelchair and whose upper limbs are partially affected and have good functionality of their arms and partial functionality of their hands, enough to accelerate, brake and turn without risk. For

these users, models with the TETRA set-up are indicated. For those with spinal injuries, the most severe level that can use models with the TETRA set-up is functional C6. It is essential to have good biceps function and firm wrist extension.

- With a total or partial disability affecting the lower limbs, who use a wheelchair and whose upper limbs are partially affected and have full functionality (arm and hand) on at least one side. For these users, models with the HEMI set-up are indicated.

Additionally, for safety reasons, the HANDBIKE BATEC ELECTRIC may only be used by people who:

- Have been instructed on how to operate it.
- Have the physical and psychological conditions required to operate it safely in all situations.
- Are able to brake the wheelchair and stop it if the device fails.

CONTRAINDICATIONS

The HANDBIKE BATEC ELECTRIC is not recommended in the following cases:

- People whose physical and/or psychological condition does not allow them to move around safely in situations with traffic.
- People who weigh over 110 kg.
- For the quad models, quadriplegic people who don't have at least the arm mobility described in this user manual.
- For the standard models, people who don't have enough hand strength to stop the device using the brake levers.
- For the hemi models, people who don't have enough strength in their working hand to stop the device using the brake levers.
- People with visual impairment that impedes them from legally driving a vehicle.

If you are unsure, please check with a professional.

SPEED

- The maximum allowed speed for the AIR model is 15 km/h, in pedestrian areas 6 km/h is recommended.
- The maximum allowed speed for the PLUS range is 20 km/h, in pedestrian areas 6 km/h is recommended.
- The maximum speed for the SPORT range is 30 km/h, in pedestrian areas 6 km/h is recommended.

However, users are responsible for knowing and obeying the speed limits on different types of roads regulated by the country, region, or locality. They should also be aware of any insurance requirements, driver's license requirements, and the need for protective and visibility equipment.

2. SAFETY MEASURES AND GENERAL PRECAUTIONS

2.1. BEFORE USING YOUR BATEC HANDBIKE FOR THE FIRST TIME

IMPORTANT!

Before using your HANDBIKE BATEC ELECTRIC for the first time, read this user manual carefully.

- Do not exceed the maximum weight of 110 kg. This includes the user's weight and any other object being carried.
- The HANDBIKE BATEC ELECTRIC isn't compatible with all wheelchairs on the market. Before installing it, check the current BATEC Compatibility Table Link <https://support.batec-mobility.com/es/support/home>.
- It is important to start out slowly the first time you use your HANDBIKE BATEC ELECTRIC in order to learn the sensations you get and become familiar with it. Each new manoeuvre is a potential risk. Be particularly careful with turns and speed.
- Choose a calm place to use your HANDBIKE BATEC ELECTRIC the first time.
- Have someone help you become familiar with your HANDBIKE BATEC ELECTRIC. Each new manoeuvre can put you at risk of injury. Never use your HANDBIKE BATEC ELECTRIC alone until you have gained the necessary skill level. When someone is helping you, it is important for them to have read the user manual beforehand as well.
- Safe use of the HANDBIKE BATEC ELECTRIC AIR and PLUS MODELS requires basic training and experience.
- On the other hand, the safe use of the HANDBIKE BATEC ELECTRIC SPORT models requires more specific training and experience than the other models.
- The HANDBIKE BATEC ELECTRIC is a medical device classified as a technical aid for mobility and we recommend consulting your local regulations on the circulation of such vehicles. BATEC MOBILITY S.L. can provide you with the Declaration of Conformity as a medical device or certification as a technical aid, if required.
- Check local laws regarding add-on handbikes and wheelchairs, as they may vary depending on your location.
- Always do a safety check on your HANDBIKE BATEC ELECTRIC before each use. To do this, check:
 - Conditions of the fixed or removable axles on the rear wheels of your wheelchair.
 - The Velcro on the seat and backrest of your wheelchair.

- Tyre pressure on your wheelchair. If it is too low, there is a greater risk of tipping over and greater wear (recommended pressure for wheelchair: see side of tyre).
- HANDBIKE BATEC ELECTRIC tyre pressure (recommended pressure: 2 bar).
- Conditions of the HANDBIKE BATEC ELECTRIC frame.
- All parts, in particular any screws, are properly tightened.
- Brakes are working properly.
- Make sure that all electrical connections are in good conditions and properly connected.
- Make sure the battery is connected.
- The HANDBIKE BATEC ELECTRIC attaches and detaches from the wheelchair properly.

2.2. DURING USE

- The HANDBIKE BATEC ELECTRIC must only be used to transport one person at a time.
- The HANDBIKE BATEC ELECTRIC is only recommended for people with the physical and mental state necessary to travel in situations with vehicle traffic.
- Avoid unnecessary risks like excess speed, fast turns and leaning to the side. Also avoid travelling over objects that may cause you to brake sharply or tip over.
- Any damage resulting from jumping while using the product is excluded from the warranty.
- The safest way to drive the HANDBIKE BATEC ELECTRIC is on level surfaces without obstacles. Use ramps whenever possible. Always be extremely careful when crossing the road.
- Use reflectors on your HANDBIKE BATEC ELECTRIC and wheelchair or on your clothes to make you more visible to drivers. In situations with reduced visibility, always use lights.
- There is a risk of skidding on wet surfaces, gravel or irregular terrain, reduce your speed and use extra caution. Do not use the HANDBIKE BATEC ELECTRIC on mud or ice.
- Make sure your feet are properly positioned on the footrest of the wheelchair and that the seatbelt is properly adjusted.
- When you use your HANDBIKE BATEC ELECTRIC, control your speed and only ride as fast as appropriate for the given ground conditions. By maintaining a prudent speed, you will avoid unnecessary risks. Don't ride the HANDBIKE BATEC ELECTRIC faster than you would travel in a wheelchair alone.
- The usual maximum speed on sidewalks and pedestrian areas is 6 km/h.

- When you need to turn, remember to slow down and be extra careful. If you take a turn too quickly, you may lose control of the HANDBIKE BATEC ELECTRIC and run the risk of tipping over or crashing. Lean your body into the curve for better stability and to reduce the risk of tipping over.
- While driving and braking, hold the grips firmly with both hands in order to avoid any accidents caused by steering error.

CAREFUL!

- Going up or down kerbs isn't recommended, as there is a risk of tipping over. When you can't avoid a kerb, we recommend you:
 - Get help whenever possible.
 - Take the kerb perpendicularly, so that both wheels of your wheelchair go over the kerb at the same time. Going over a kerb in any other way significantly increases the risk of tipping over.
 - Avoid banging the HANDBIKE BATEC ELECTRIC when going up or down over a kerb, as this could damage both the HANDBIKE BATEC ELECTRIC and your wheelchair.
 - The HANDBIKE BATEC ELECTRIC and the wheelchair must roll gently over the kerb on all wheels.
 - See also figures 61 and 62.
-
- When travelling up or down a hill, we recommend you:
 - Get help whenever possible while on a hill.
 - Stay a straight course.
 - When going up steep hills, lean forward to put more weight on the front wheel and improve traction.
 - Avoid unnecessary turns and never turn sharply, as this puts you at risk of tipping over.
 - When going downhill, go slowly, braking continuously, and avoid braking sharply.
 - Don't use the brakes on your wheelchair; this could cause them to lock up and you to lose control.
 - The maximum recommended slope is 10°, as steeper inclines could make the parking brake less effective. Also take into account that the safe braking distance on hills can be much longer than on flat terrain.
 - When travelling long distances we recommend you carry an extra battery.
 - Make sure the grip on your wheels is proper conditions.
 - Make sure you are familiar with braking distances at various speeds.

- Only use accessories approved by BATEC MOBILITY.
<https://batec-mobility.com/en/batteries-and-accessories/>
- NEVER use your HANDBIKE BATEC ELECTRIC with a child on your lap, as there is a high risk
- The following manoeuvres are not allowed:
 - Riding up or down stairs or staircases.
 - Riding at an angle on the pavement.
 - 180° turns on steep hills.
 - Sharp movements.
 - Sharp braking on a curve.

CAREFUL!

There is a risk that fingers may get smashed in several parts of the product (openings and spaces between fixed accessible parts) during intended use. (See fig. 1, 2)

CAREFUL!

There is a risk that feet may get smashed (openings and spaces between fixed accessible parts) during intended use. (See fig. 3)

2.3. GENERAL PRECAUTIONS

CAREFUL!

- Remember to always turn off the BATEC ELECTRIC HANDBIKE before attaching or detaching it to the wheelchair.
- Do not use the BATEC ELECTRIC HANDBIKE if it behaves abnormally or erratically.
- Do not use the BATEC ELECTRIC HANDBIKE with dead batteries.
- Take your HANDBIKE BATEC ELECTRIC to a distributor or technical service provider authorised by BATEC MOBILITY regularly for servicing.
- Avoid contact or proximity to fire, particularly lit cigarettes.
- If the HANDBIKE BATEC ELECTRIC is exposed to direct sunlight or low temperatures for an extended period of time, the HANDBIKE BATEC ELECTRIC pieces may overheat (>41°C) or become too cold (<0 °C).
- Only use our charger to charge the battery.

- If you don't use your HANDBIKE BATEC ELECTRIC for an extended period of time, make sure the battery is fully charged and recharge it every 2 months to avoid damage caused by "deep discharge".
- Keep the HANDBIKE BATEC ELECTRIC out of the reach of children when not under adult supervision.
- Do not open the battery, as there is a risk of explosion and/or fire.
- The repetitive rotational movement when operating the throttle grip may cause a wrist injury known as intersection syndrome. It is recommended to rest the wrist and avoid intensive use of the throttle grip.
- Do not tamper with the connector located at the bottom of the battery holder. Removing the protective cover could cause the controller to malfunction if the connector is exposed to adverse conditions.

2.4. LIFESPAN

The BATEC ELECTRIC HANDBIKE has a lifespan of at least 5 years, as long as it is used according to the indicated operating conditions and receives proper maintenance and revisions, as indicated.

The life of the battery is between 200 and 500 charging cycles, and between 1 and 3 years, depending on several factors, including temperature, frequency of use, consistency of use and other individual issues (user weight, driving style, terrain, etc.).

Disposing of the product

The BATEC ELECTRIC HANDBIKE is an electric/electronic device, so it must be disposed of properly. Check with your national distributor or a local waste-management centre if you have any questions.

3. MODELS: CHARACTERISTICS AND PARTS

i INFORMATION!

The Batescs shown and described in this manual may not be exactly the same as the model you own. However, all of the instructions provided here apply. The manufacturer reserves the right to modify the weights, measurements and other technical details in this manual without prior notice. All figures, measurements and capacities in this manual are approximate and not technical specifications. The manufacturer reserves the right to modify this manual without prior notice. The most current version of the manual is always available online: <https://support.batec-mobility.com/en/support/>

There are several models of HANDBIKE BATEC ELECTRIC: AIR / MINI 2 PLUS, ELECTRIC 2 PLUS and SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT. All of the HANDBIKE BATEC ELECTRIC models have the same anchoring system and the same type of structure, only wheel size, top speed and features:

- AIR: Powered by an electric motor installed on the 10" front wheel. (See fig. 4a)
- MINI 2 PLUS: powered by an electric motor on the 12" front wheel. (See fig. 4b)
- ELECTRIC 2 PLUS, RAPID 2 PLUS and RAPID SPORT: powered by an electric motor on the 18" front wheel. (See fig. 4c)
- SCRAMBLER 2 PLUS and SCRAMBLER SPORT : powered by an electric motor on the 19" front wheel (all terrain). (See fig. 4d)

The various models of HANDBIKE BATEC ELECTRIC are designed for users with full hand function in their standard configuration. For users without full hand function, we have specially adapted controls, sold under the name QUAD. For users who only have functional difficulties in one of their upper limbs, we offer controls specially adapted to this situation, sold under the name HEMI. More information in the box below. The HANDBIKE BATEC ELECTRIC models (AIR, MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT) are sold with different speed limits depending on the country, region or client. (See fig. 60)

i INFORMATION!

SKILLS NEEDED TO USE THE BATEC HANDBIKES: HANDBIKES BATEC are only recommended for people with the physical and psychological conditions to circulate in traffic. The standard HANDBIKE BATEC ELECTRIC models require full use of hands and arms.

To use the models with QUAD controls, users must have good arm function and partial hand function, enough to accelerate, brake and turn without danger. This is why BATEC MOBILITY always offers users the option to test the product before purchasing.

The most serious spinal cord injury compatible with use of the HANDBIKE BATEC ELECTRIC QUAD is functional C6 (good biceps function and wrist motion is essential). To use versions equipped with HEMI controls, users must have full functionality (hand and arm) in at least one of their upper limbs. (See fig. 5)

3.1. HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT

3.1.1. CHARACTERISTICS

AIR

- Material: Aluminium 7005 T6 / Aluminium 6061 T6
- Motor: Brushless 350 W 350 RPM / Peak power: 470 W / Reverse gear
- Top speed: Depending on the version. (See fig. 60)
- Power use: 10.8 Wh/km (power use calculated under the following conditions: 75 kg user, extra-light wheelchair, flat terrain, no wind and steady speed of 10 km/h. May vary $\pm 20\%$ depending on the battery).
- 2A charger
- The flame resistance of the materials and polymeric components of the HANDBIKE BATEC ELECTRIC power and steering system is V-0, according to UL 94 standards
- Tyre: 10" urban
- Rim: 10" aluminium
- Brakes: Dual mechanical brake callipers with a 160 mm disc and aluminium parking brake levers
- Headset: 1.5" black
- Handlebar: Black matte 6061-T6 aluminium
- HANDBIKE BATEC ELECTRIC EASY-FIT regulation system (pat. pend.)
- HANDBIKE BATEC ELECTRIC EASY-FIX anchor system (pat. pend.)
- HANDBIKE BATEC ELECTRIC SAFE-FEET stand frame (pat. pend)
- 30 lux front LED light, runs off main battery

- Dual LED tail lights on the stand frame
- TFT COLOT ALL IN ONE display screen (battery gauge, 3-speed forward selector and 1 reverse, 3 assistance modes, cycle computer, front and tail lights and horn, controlled by 3 buttons on the handlebar)
- Weight: 11 kg (without battery)
- Dimensions: (See fig. 5A, 5B and 5C)

MINI 2 PLUS

- Material: Aluminium 7005 T6/ Aluminium 6061 T6
- Motor: Brushless 350 W 350 RPM / Peak power: 600 W / Reverse gear
- Top speed: depending on the model. (See fig. 60)
- Power use: 13 Wh/km (power use calculated under the following conditions: 75-kg user, extra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. May vary ± 20 % depending on the battery).
- Super-fast 4A charger
- The BATEC power and steering system has a V-0 flammability rating, according to UL94 standards.
- Tyre: 12" ultra-grip
- Rim: 12" aluminium
- Brakes: 140 mm mechanical disc brakes with aluminium levers and parking brake.
- Headset: 1.1/8 black
- Handlebar: 7005-T6 / 6061-T6 matte black aluminium
- BATEC EASY-FIT regulation system (Pat. Pend.)
- BATEC EASY-FIX anchor system (Pat. Pend.)
- BATEC SAFE-FEET stand frame (Pat. Pend.)
- 700 lumen LED front light connected to the main battery
- Dual LED tail light on the stand frame
- ALL-IN-ONE display panel (battery indicator, 3-speed selector and rear gear, cycle computer, front and rear lights, and horn, all with 4 buttons on the handlebar)
- Weight: 13.9 kg (without battery)
- Dimensions: (See fig. 6, 7)

ELECTRIC 2 PLUS

- Material: 7005-T6 / 6061-T6 aluminium / reinforced removable dropout to change the tyre

- Motor: 250 RPM, 900 W brushless motor / Nominal power: 900 W / Reverse gear
- Top speed: depending on the model. (See fig. 60)
- Power use: 13 Wh/km (power use calculated under the following conditions: 75-kg user, extra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. Life may vary ± 20 % depending on the battery).
- Fast 4A charger
- Tyre: Exclusive ultra-grip 18"
- Rim: 18" double-wall aluminium
- Brakes: 180 mm mechanical disc brakes. Aluminium levers with parking brake and throttle cut-off.
- Headset: 1.1/8 black
- Handlebar: 6061-T6 black matte aluminium
- HANDBIKE BATEC ELECTRIC EASY-FIT regulation system
- HANDBIKE BATEC ELECTRIC EASY-FIX anchor system
- HANDBIKE BATEC ELECTRIC SAFE-FEET stand frame
- 700 lumen front LED light, runs off main battery
- Dual LED tail light on the stand frame
- ALL IN ONE LCD display (battery level indicator, 3 forward and 1 reverse speed selector, cruise control mode indicator and anti-slip mode, cycle computer, front and rear lights, and horn; operated with 6 buttons on the handlebars)
- Removable weights
- Weight: 14.9 kg (without battery/weights)
- Dimensions: (See fig. 8)

RAPID 2 PLUS / RAPID SPORT

- Colour: graphite with red trim
- Top speed: depending on the model. (See fig. 60)
- Power use: 13 Wh/km (power use calculated with: 75-kg user, ultra-light wheelchair, flat terrain, no wind and a steady speed of 15 km/h. May vary ± 20 % depending on the battery).
- Motor: 300-RPM brushless / Nominal power: 1440 W / Peak power: 1620 W
- Reverse gear
- Batteries with connector only suitable for use with the BATEC RAPID 2 PLUS and BATEC SCRAMBLER 2 PLUS 48 V
- Fast 3A charger

- Chassis: Black aluminium 7005 T6 headset 1.1/8 / Removable dropout to change the tyre
- 6061-T6 matte black aluminium handlebar
- Mechanical brakes / Aluminium levers with parking brake and throttle cut-off. 180 mm floating rotor disc
- Tyre: Exclusive ultra-grip 18"
- Rim: 18" black double-wall aluminium
- 700-lm front light with running lamp. Dual LED tail light on the stand frame
- ALL IN ONE LCD display (battery level indicator, 3 forward and 1 reverse speed selector, cruise control mode indicator and anti-slip mode, cycle computer, front and rear lights, and horn; operated with 6 buttons on the handlebars)
- USB port
- Horn
- BATEC ACOPLA support
- Removable Scrambler QR weights
- Weight 15.9 kg (without battery or weights)
- Dimensions: (See fig. 9)

SCRAMBLER 2 PLUS / SCRAMBLER SPORT

- Material: 7005-T6 / 6061-T6 aluminium / reinforced removable dropout to change the tyre
- Motor: 250 RPM, 900 W brushless motor / Nominal power: 1440 W / Peak power: 1620 W / Reverse gear
- Top speed: depending on the model. (See fig. 60)
- 48V batteries with connector only suitable for use with the BATEC RAPID 2 PLUS and BATEC SCRAMBLER 2 PLUS
- Power use: 13 Wh/km (power use calculated under the following conditions: 75-kg user, extra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. Life may vary ± 20 % depending on the battery).
- Fast 3A charger
- Tyre: 19" ultra-grip
- Rim: 19" double-wall aluminium
- Brakes: mechanical brakes / Aluminium levers with parking brake and throttle cut-off / 203-mm rotor disc
- Headset: 1.1/8 FSA Orbit Mx black
- Handlebar: 6061-T6 black matte aluminium

- HANDBIKE BATEC ELECTRIC EASY-FIT regulation system
- HANDBIKE BATEC ELECTRIC EASY-FIX anchor system
- Stand frame with suspension
- New dual 1400-lumen front lights
- Dual LED tail light on the stand frame
- ALL IN ONE LCD display (battery level indicator, 3 forward and 1 reverse speed selector, cruise control mode indicator and anti-slip mode, cycle computer, front and rear lights, and horn; operated with 6 buttons on the handlebars)
- Removable Scrambler QR weights
- Weight: 16.1 kg (without battery/weights)
- Dimensions: (See fig. 10)

3.1.2. PARTS OF THE HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT (See fig. 11a, 11b, 11c, 11d, 11e, 11f, 11g and 11A, 11B, 11C, 11D, 11E, 11F, 11G for Air)

3.1.3. PARTS OF THE TETRA CONTROL FOR THE HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT (See fig. 11f, 11g, 12, and 12A for Air)

3.1.4. PARTS OF THE HEMI CONTROL FOR THE HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT (See fig. 11f, 11g, 13 and 13A for Air)

A left-handed version is also available with symmetrical controls regarding the one displayed.

4. OPERATING THE HANDBIKE BATEC ELECTRIC / USER GUIDE

4.1. DELIVERY

Your new HANDBIKE BATEC ELECTRIC will be delivered in a cardboard box. In order to prevent damage during transport, parts that must be assembled later, like the battery, are packaged separately inside the cardboard box.

4.2. UNPACKING

- Check that the outside of the cardboard box hasn't been damaged in any way during transport.
- Remove all of the packing materials.
- Remove the parts of the HANDBIKE BATEC ELECTRIC (including the battery and charger) from the cardboard box and check that they are all there.
- Check that all the parts are in proper conditions, without any sort of damage, scratches, tears, dents, irregularities or other defects.
- If you find any problems, or if any pieces are missing, notify the authorised distributor immediately.

4.3. ELEVATION, TRANSPORT AND STORAGE

The HANDBIKE BATEC ELECTRIC has been designed so that users can easily remove and detach the battery on their own.

- When taking the HANDBIKE BATEC ELECTRIC on any means of transport (car, bus, train, plane, etc.), we recommend removing the battery before your journey: turn the battery key to OPEN, remove the key and then remove the battery.
- When lifting the HANDBIKE BATEC ELECTRIC, we recommend removing the weight as indicated in section 3.10 and the battery as per section 4.9.1.
- When transporting the HANDBIKE BATEC ELECTRIC in the boot of a vehicle, we recommend leaving it attached to a wheelchair whenever possible.

The BATEC HANDBIKE is not designed to be used as a means of transport inside a vehicle.

The HANDBIKE BATEC ELECTRIC ELECTRIC should be stored in a cool, dry place, protected from direct sunlight and rain. It should not be kept out of use for more than two months.

4.4. ATTACHING THE HANDBIKE BATEC ELECTRIC:

CAUTION!

Before anchoring the Batec Electric Air handbike, make sure that the clamp levers are secure. When tightening the clamp levers, when they begin to resist, rotate 1/8 of a turn of retightening to ensure that the travel is locked.

- Step 1: Get in front of the HANDBIKE BATEC ELECTRIC, with the handlebar centred in front of the wheelchair, with the brake on. (See fig. 14a, 14b)
- Step 2: Grab the HANDBIKE BATEC ELECTRIC by the handlebar and pull it towards you, rolling the frame stand until the detaching ramp on the anchoring bar touches the front axle of the bar on the wheelchair. (See fig. 15)
- Step 3: Tip the handlebar, lifting the front wheel, while pulling the HANDBIKE BATEC ELECTRIC towards you. When the detaching ramp has passed the front axle of the bar on the wheelchair, let go of the handlebar and put the front wheel back on the ground. (See fig. 16a, 16b)
- Step 4: Take off the brake on the wheelchair and push gently on the HANDBIKE BATEC ELECTRIC handlebar to make sure the front axle is securely in place.
- Step 5: Push the handle until it clicks, so that both the wheels on the HANDBIKE BATEC ELECTRIC stand frame and the front wheels on the wheelchair are off the ground. Your HANDBIKE BATEC ELECTRIC is now attached. (See fig. 17)

4.5. DETACHING THE HANDBIKE BATEC ELECTRIC

- Step 1: Make sure your HANDBIKE BATEC ELECTRIC is switched off. (See fig. 18 or 18A for AIR)
- Step 2: Push the anchor lever forward with the palm of your hand until you feel it come free. Don't let go of the lever and set the HANDBIKE BATEC ELECTRIC on the ground. Advice: lean your torso back to take weight off the lever. This will make it much easier to detach the HANDBIKE BATEC ELECTRIC. (See fig. 19 or 19A for AIR)
- Step 3: Tilt the handlebar, lifting the front wheel, while you push the HANDBIKE BATEC ELECTRIC away from you. (See fig. 20)
- Step 4: When the detaching ramp has gone over the front axle of the wheelchair bar, let go of the handlebar and set the front wheel down again. Your HANDBIKE BATEC ELECTRIC is now detached. (See fig. 21)

.....

 **ADVISE!**

Practice attaching and detaching your HANDBIKE BATEC ELECTRIC calmly and you'll see that, in just a few weeks of use, it will become automatic

.....

4.6. FOLDING AND UNFOLDING THE HANDBIKE BATEC ELECTRIC AIR AND MINI 2 PLUS

The MINI 2 PLUS and AIR models are the only ones that can be folded for easier transport when inactive.

4.6.1 FOLDING THE HANDBIKE BATEC ELECTRIC AIR ANCHORING

- Step 1. Loosen the union levers connecting the anchoring bar to the chassis by one turn. (See fig. 21A)
- Step 2. Slide the anchoring bar to the highest point of the guide. (See fig. 21B)
- Step 3. Rotate the anchoring bar until its other end makes contact with the fender/foot protector. (See fig. 21C)
- Step 4. Lower the handlebar to its minimum height and rotate it 90 degrees to reduce the HANDBIKE BATEC ELECTRIC AIR to its most compact form. (See fig. 21D)
- Step 5. Re-tighten the lower clamping lever to restrict the rotation of the anchoring bar and prevent accidents.

4.6.2. UNFOLDING THE HANDBIKE BATEC ELECTRIC AIR ANCHORING

.....

 **CAUTION!**

During unfolding process, ensure that the cables do not get trapped with other components of the device.”

.....

- Step 1. Loosen the screw that restricts the rotation of the anchoring bar.
- Step 2. Once the anchoring bar can rotate freely, slightly raise it to its limit. (See fig. 21E)
- Step 3. Rotate the anchoring bar until it hits the stop, which indicates the bar is horizontal. (See fig. 21F)
- Step 4. Lower the anchoring bar, ensuring it follows the movement guide until it makes contact with the stop that secures the connection height of the anchoring bar with the wheelchair. (See fig. 21G)

- Step 5. Tighten the CLAMPING LEVERS of the anchoring bar against the chassis, ensuring the bar is fully secured. When tightening, once the CLAMPING LEVERS start to resist, turn an additional 1/8 of a turn to ensure anchoring. (See fig. 21H)
- Step 6. By pressing the CLAMPING LEVERS inward and rotating them on their axle, you can adjust them to your preferred position.

4.6.3. FOLDING THE HANDLEBAR OF THE HANDBIKE BATEC ELECTRIC MINI 2 PLUS

- Step 1. Turn the handlebar 45° to the left.
- Step 2. Rotate the safety latch clockwise to unlock the release lever. (See fig. 21I)
- Step 3. Engage the release lever to free the locking shaft of the folding hinge. (See fig. 21J)
- Step 4. Fold the handlebar to the right. (See fig. 21K)
- Step 5. Secure the handlebar against the stand frame using the elastic cord. (See fig. 21L)

4.6.4. UNFOLDING THE HANDLEBAR OF THE HANDBIKE BATEC ELECTRIC MINI 2 PLUS

- Step 1. Release the elastic cord.
- Step 2. Raise the handlebar.
- Step 3. Engage the locking lever to secure the folding hinge of the handlebar.
- Step 4. Rotate the safety latch to lock the locking lever in place.

4.7. USING THE HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT CONTROLS

You can find out the meaning of all the symbols used on the controls of the HANDBIKE BATEC ELECTRIC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT in Figure 58, and for the AIR in Figure 58A, with their respective positions in Figures 58B, 58C, 58D, 58E and 58F.

4.7.1. TURN ON

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key on the right side of the battery clockwise to the ON position. It will do a system check: don't start using your HANDBIKE BATEC ELECTRIC until it has finished with the system check and all the lights are on. (See fig. 22)

→ IMPORTANT!

For the safety check to be completed properly, the wheel of the HANDBIKE BATEC ELECTRIC must be stationary and the user must press one of the BRAKE LEVERS at least once. (See fig. 22B and 22C)

BATEC AIR

- Before switching on the battery, ensure it is properly connected and secured with the key.
- Press the "I" button on the battery to turn on the HANDBIKE BATEC ELECTRIC AIR.
- Then, press the "M" button on the display panel to switch on the TFT screen. (See fig. 22A)

→ IMPORTANT!

For the safety check to be completed properly, the wheel of the HANDBIKE BATEC ELECTRIC must be stationary and the user must press one of the BRAKE LEVERS at least once. (See fig. 22D and 22E)

4.7.2. TURN OFF

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key on the battery anti-clockwise to the OFF position. If you leave it in the OPEN position, your HANDBIKE BATEC ELECTRIC will be off and the battery will be ready to remove. (See fig. 22)

BATEC AIR

- Press the "O" button on the battery to turn off the BATEC HANDBIKE ELECTRIC.
- Do not remove the battery if the switch is in the "I" position. (See fig. 22A)

4.7.3. SPEED UP

The throttle is on the right grip on the handlebar. It works just like a motorbike, meaning you can regulate your speed by turning the throttle on the grip. (See fig. 23)

4.7.4. BRAKE

The HANDBIKE BATEC ELECTRIC has two brake levers on the handlebar that activate the disc brakes. Squeeze the lever to activate the brakes. Advice: Use one hand to accelerate and the other to brake. Advice: Always keep two fingers on the

brake lever in order to respond quickly. Pressing on the brake lever and the red lever at the same time blocks the brake lever, acting as a parking brake. (See fig. 24 or 24A for AIR). The MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT models feature an electronic brake that helps maintain a constant speed. When going downhill, the electronic brake prevents the vehicle from exceeding the maximum speed.

4.7.5. FRONT AND BACK GEARS

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

- In front gear, there are 3 speed limits (position 1, position 2 and position 3), each with a light indicator on the display panel.
- In back gear, there is only one speed and one light indicator on the display panel.
- When the HANDBIKE BATEC ELECTRIC is switched on, it will always be in front gear, position 1.
- In front gear, press button F to choose the desired speed limit (position 1, position 2 or position 3).
- Press the R button to switch to back gear. This mode will only be activated if the HANDBIKE BATEC ELECTRIC is stopped and the throttle is not engaged. Press the F button to switch back to front gear. (See fig. 25)

BATEC AIR

- For forward motion, there are 3 available speed limits (position 1, position 2, position 3/4), each represented by a different image on the screen.
- For reverse motion, there is only one speed limit, shown as "R" on the screen.
- When the HANDBIKE BATEC ELECTRIC is switched on, it will always be in forward motion, position 1.
- Pressing the buttons ▲ or ▼ allows you to select the desired speed limit (position 1, position 2, position 3, position 4 or reverse "R").
- Pressing the button ▼ activates reverse mode. This mode can only be engaged if the HANDBIKE BATEC ELECTRIC is nearly stationary and the throttle is disengaged. Pressing the button ▲ returns the system to forward mode if the HANDBIKE BATEC ELECTRIC is almost stationary and the throttle is disengaged. (See fig. 25A)
- There are three operating modes: ECO, NORMAL and POWER
 - Positions "R", 1 and 2 are ECO
 - Position 3 is NORMAL
 - Position 4 is POWER

4.7.6. FREEWHEEL MODE

Turn off the device using the key (BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS) or turn off the battery switch (BATEC AIR)

- To move forwards or backwards, use the rims on the wheelchair to move the HANDBIKE BATEC ELECTRIC and wheelchair in the desired direction.
- To turn, turn the HANDBIKE BATEC ELECTRIC with one hand while propelling the wheelchair with the other.
- To stop, use the HANDBIKE BATEC ELECTRIC brakes or the rim on the wheelchair.

4.7.7. LIGHTS

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

There is a button on the HANDBIKE BATEC ELECTRIC handlebar to switch on the front and rear lights. When you push the button, the front light will go on in normal mode and the light indicator on the display panel will turn green. The rear lights have a flashing function when braking while the lights are on. When you push the front light button again, it changes to high-beam mode, and so on. To turn off the lights, hold the button down for 2 seconds. The front light is below the battery and the two rear lights, on either side of the stand frame. (See fig. 26a, 26b, 26c)

BATEC AIR

To turn on the front and tail lights of the HANDBIKE BATEC ELECTRIC, press and hold the button  for 1 second. When the lights are on, the brightness of the LCD screen will dim to prevent visual discomfort during night-time use. To turn off the lights, press and hold the button  for 1 second, and the LCD screen brightness will return to its normal level. The front light is below the battery and the two tail lights, on either side of the stand frame. (See fig. 26A, 26B, 26c)

4.8. USING THE HANDBIKE BATEC ELECTRIC QUAD AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT CONTROLS

You can find out the meaning of all the symbols used on the controls of the HANDBIKE BATEC ELECTRIC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT in Figure 58, and for the AIR in Figure 58A, with their respective positions in Figures 58B, 58C, 58D, 58E and 58F.

4.8.1. TURN ON

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key on the right side of the battery clockwise to the ON position. It will do a system check: don't start using your HANDBIKE BATEC ELECTRIC until it has finished with the system check and all the lights are on. (See fig. 27)

IMPORTANT!

For the safety check to be completed properly, the wheel of the HANDBIKE BATEC ELECTRIC must be stationary and the user must apply the brake lever on the throttle side at least once. (See fig. 22B and 22C)

BATEC AIR

- Before switching on the battery, ensure it is properly connected and secured with the key.
- Press the "I" button on the battery to turn on the HANDBIKE BATEC ELECTRIC AIR.
- Then, press the "M" button on the display panel to switch on the TFT screen. (See fig. 27A)

IMPORTANT!

For the safety check to be completed properly, the wheel of the HANDBIKE BATEC ELECTRIC must be stationary and the user must press one of the BRAKE LEVERS at least once. (See fig. 22D and 22E)

4.8.2. TURN OFF

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key on the battery anti-clockwise to the OFF position. If you leave it in the OPEN position, your HANDBIKE BATEC ELECTRIC will be off and the battery will be ready to remove. (See fig. 27)

BATEC AIR

- Press the "O" button on the battery to turn off the BATEC HANDBIKE ELECTRIC.
- Do not remove the battery if the switch is in the "I" position. (See fig. 27A)

4.8.3. SPEED UP

Pull the HANDBIKE BATEC ELECTRIC QUAD MINI handlebar towards you to activate the throttle. (See fig. 28)

4.8.4. BRAKE

Push the HANDBIKE BATEC ELECTRIC QUAD MINI handlebar away from you to activate the brake. (See fig. 29)

4.8.5. FRONT AND BACK GEARS

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORTS

- In front gear, there are 3 speed limits (position 1, position 2 and position 3), each with a light indicator on the display panel.
- In back gear, there is only one speed and one light indicator on the display panel.
- When the HANDBIKE BATEC ELECTRIC is switched on, it will always be in front gear, position 1.
- In front gear, press button F to choose the desired speed limit (position 1, position 2 or position 3).
- Press the R button to switch to back gear. This mode will only be activated if the HANDBIKE BATEC ELECTRIC is sopped and the throttle is not engaged. Press the F button to switch back to front gear. (See fig. 30)

BATEC AIR

- For forward motion, there are 3 available speed limits (position 1, position 2, position 3/4), each represented by a different image on the screen.
- For reverse motion, there is only one speed limit and an indicator light on the display panel.
- When the HANDBIKE BATEC ELECTRIC is switched on, it will always be in forward motion, position 1.
- Pressing the buttons  or  allows you to select the desired speed limit (position 1, position 2, position 3, position 4 or reverse "R").
- Pressing the button  activates reverse mode. This mode can only be engaged if the HANDBIKE BATEC ELECTRIC is nearly stationary and the throttle is disengaged. Pressing the button  returns the system to forward mode if the HANDBIKE BATEC ELECTRIC is almost stationary and the throttle is disengaged. (See fig. 30A)
- There are three operating modes: ECO, NORMAL and POWER
 - Positions "R", 1 and 2 are ECO
 - Position 3 is NORMAL
 - Position 4 is POWER

4.8.6. FREEWHEEL MODE

To use the HANDBIKE BATEC ELECTRIC in freewheel mode:

- Turn off the device using the key (BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT) or turn off the battery switch (BATEC AIR).- Switch the device off with the key.
- To move forwards or backwards, use the rims on the wheelchair to move the HANDBIKE BATEC ELECTRIC and wheelchair in the desired direction.
- To turn, turn the HANDBIKE BATEC ELECTRIC with one hand while propelling the wheelchair with the other.
- To stop, use the HANDBIKE BATEC ELECTRIC brakes or the rim on the wheelchair.

4.8.7. LIGHTS

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

There is a button on the HANDBIKE BATEC ELECTRIC handlebar to switch on the front and rear lights. When you push the button, the front light will go on in normal mode and the light indicator on the display panel will turn green. The rear lights have a flashing function when braking while the lights are on. When you push the front light button again, it changes to high-beam mode, and so on. To turn off the lights, hold the button down for 2 seconds. The front light is below the battery and the two rear lights, on either side of the stand frame. (See fig. 31a, 31b, 31c)

BATEC AIR

To turn on the front and tail lights of the HANDBIKE BATEC ELECTRIC, press and hold the button  for 1 second. When the lights are on, the brightness of the TFT screen will dim to prevent visual discomfort during night-time use. To turn off the lights, press and hold the button  for 1 second again, and the LCD screen brightness will return to its normal level. The front light is below the battery and the two tail lights, on either side of the stand frame. (See fig. 31A, 31B, 31C)

4.9. OPERATING THE HANDBIKE BATEC ELECTRIC AIR / MINI 2 PLUS / ELECTRIC 2 PLUS (RAPID 2 PLUS) / SCRAMBLER 2 PLUS / RAPID SPORT / SCRAMBLER SPORT HEMI CONTROLS

You can find out the meaning of all the symbols used on the controls of the HANDBIKE BATEC ELECTRIC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS and SCRAMBLER 2 PLUS in Figure 58, and for the AIR in Figure 58A, with their respective positions in Figures 58B, 58C, 58D, 58E and 58F.

4.9.1. TURN ON

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key on the right side of the battery clockwise to the ON position. It will do a system check: don't start using your HANDBIKE BATEC ELECTRIC until it has finished with the system check and all the lights are on. (See fig. 32)

IMPORTANT!

For the safety check to be completed properly, the wheel of the HANDBIKE BATEC ELECTRIC must be stationary and the user must press one of the BRAKE LEVERS at least once. (See fig. 22B and 22C)

BATEC AIR

Press the "I" button on the battery to turn on the HANDBIKE BATEC ELECTRIC. Before switching it on, ensure the battery is properly connected and secured with the key. Then, press the "M" button to power on the TFT screen. (See fig. 32A)

IMPORTANT!

For the safety check to be completed properly, the wheel of the HANDBIKE BATEC ELECTRIC must be stationary and the user must press one of the BRAKE LEVERS at least once. (See fig. 22D and 22E)

4.9.2. TURN OFF

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key on the battery anti-clockwise to the OFF position. If you leave it in the OPEN position, your HANDBIKE BATEC ELECTRIC will be off and the battery will be ready to remove. (See fig. 32)

BATEC AIR

- Press the "O" button on the battery to turn off the BATEC HANDBIKE ELECTRIC.
- Do not remove the battery if the switch is in the "I" position. (See fig. 32A)

4.9.3. SPEED UP

The throttle is on the right or left grip on the handlebar. (See fig. 33 or 33A for AIR)

4.9.4. BRAKE

Pressing on the brake lever and the red lever at the same time blocks the brake lever, acting as a parking brake. (See fig. 34 or 34A for AIR). The MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT models feature an electronic brake that helps maintain a constant speed. When going downhill, the electronic brake prevents the vehicle from exceeding the maximum speed.

4.9.5. FRONT AND BACK GEARS

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

- In front gear, there are 3 speed limits (position 1, position 2 and position 3), each with a light indicator on the display panel.
- In back gear, there is only one speed and one light indicator on the display panel.
- When the HANDBIKE BATEC ELECTRIC is switched on, it will always be in front gear, position 1.
- In front gear, press button F to choose the desired speed limit (position 1, position 2 or position 3).
- Press the R button to switch to back gear. This mode will only be activated if the HANDBIKE BATEC ELECTRIC is sopped and the throttle is not engaged. Press the F button to switch back to front gear. (See fig. 35)

BATEC AIR

- For forward motion, there are 3 available speed limits (position 1, position 2, position 3/4), each represented by a different image on the screen.
- For reverse motion, there is only one speed limit, shown as "R" on the screen.
- When the HANDBIKE BATEC ELECTRIC is switched on, it will always be in forward motion, position 1.
- Pressing the buttons  or  allows you to select the desired speed limit (position 1, position 2, position 3, position 4 or reverse "R").
- Pressing the button  activates reverse mode. This mode can only be engaged if the HANDBIKE BATEC ELECTRIC is nearly stationary and the throttle is disengaged. Pressing the button  returns the system to forward mode if the HANDBIKE BATEC ELECTRIC is almost stationary and the throttle is disengaged. (See fig. 30A)
- There are three operating modes: ECO, NORMAL and POWER
 - Positions "R", 1 and 2 are ECO
 - Position 3 is NORMAL
 - Position 4 is POWER

4.9.6. FREEWHEEL MODE

To use the HANDBIKE BATEC ELECTRIC in freewheel mode:

- Switch the device off with the key.
- To move forwards or backwards, use the rims on the wheelchair to move the HANDBIKE BATEC ELECTRIC and wheelchair in the desired direction.
- To turn, turn the HANDBIKE BATEC ELECTRIC with one hand while propelling the wheelchair with the other.
- To stop, use the HANDBIKE BATEC ELECTRIC brakes or the rim on the wheelchair.

4.9.7. LIGHTS

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

The front and rear light switch is located on the handlebars of the BATEC ELECTRIC HANDBIKE. When the switch is pressed, the front light will illuminate in position light mode, and the rear lights and the console indicator light will turn green. Pressing the switch again will change the front light to high beam mode, and so on. To turn off the lights, press and hold the switch for 2 seconds. The front light is located under the battery, and the two rear lights are at the end of the kickstand. (See fig. 36a, 36b, 36c)

BATEC AIR

To turn on the front and tail lights of the HANDBIKE BATEC ELECTRIC, press and hold the button **▲** for 1 second. When the lights are on, the brightness of the TFT screen will dim to prevent visual discomfort during night-time use. To turn off the lights, press and hold the button **▲** for 1 second again, and the LCD screen brightness will return to its normal level. The front light is below the battery and the two tail lights, on either side of the stand frame. (See fig. 36A, 36B, 36C)

4.10. ANTI-SLIP SYSTEM

To activate the Anti-Slip System

This will occur if the BATEC ELECTRIC HANDBIKE is switched on, regardless of whether it is moving or stationary, and with cruise control deactivated. (see Fig. 58).

- Press and hold the forward drive button (“9” in figure 58) for 2 seconds.
 - The light indicator (1) (2) (3) on the display panel will blink.
 - The Anti-Slip System is now active.

When the Anti-Slip system is activated for 2 seconds, it prevents the BATEC ELECTRIC HANDBIKE from rolling backwards on a slope, and also greatly reduces wheel slippage when starting.

Deactivate the Anti-Slip System

This will occur if the BATEC ELECTRIC HANDBIKE is switched on, regardless of whether it is moving or stationary.

- Press and hold the forward drive button (“9” in figure 58) for 2 seconds
 - The light indicator (1) (2) (3) on the display panel lights up.

With the Anti-Slip system deactivated, the BATEC ELECTRIC HANDBIKE returns to its basic operation.

The Anti-Slip System is standard on the MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, and SCRAMBLER 2 PLUS models. For other models, please inquire.

4.11. USING THE CRUISE CONTROL SYSTEM

Switching on the cruise control system

This can only be done when the BATEC ELECTRIC HANDBIKE is in motion with the accelerator activated (see fig. 58).

- Hold down the **▲** button or the **▼** button (“7” or “8” on figure 58) for 2 seconds.
 - Lights (1), (2) and (3) will flash alternately on the display panel.
 - The cruise control is on.
 - The HANDBIKE will maintain the set speed shown on the screen.
- After switching on the cruise control, there are two options to adjust the set speed:
 - Press the **▲** button (“7” on figure 58) to increase the set speed.
 - Press the **▼** button (“8” on figure 58) to decrease the set speed.
 - Hold down the accelerator to increase the set speed, although when the accelerator is released the HANDBIKE will return to the previous set speed.

Switching off the cruise control system

- To switch off the cruise control system, there are three options:
 - Press the F button (“9” on figure 58).
 - Press the R button (“10” on figure 58).
 - Squeeze either of the brake levers.
- The cruise control system is switched off if:
 - Lights (1) (2) (3) on the display panel light up.
 - The HANDBIKE tends to slow down gradually until it comes to a stop if the accelerator isn’t activated.

With the cruise control off, the BATEC ELECTRIC HANDBIKE will return to its standard operating mode.

The cruise control system comes standard on the MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS and SCRAMBLER 2 PLUS models. For other models, please check.

4.12. BATTERY

BATEC MOBILITY offers several 36 V and 48 V lithium-ion batteries, all without memory effect and with a cable-free connection system.

BATEC MOBILITY offers 36 V BATTERIES THAT COMPLY WITH IATA-DGR REGULATIONS for air transport of batteries in mobility aid devices. (See fig. 37, and for AIR 37A), for the BATEC AIR, MINI 2 PLUS and ELECTRIC 2 PLUS models.

Lithium ion batteries have many advantages over other batteries, mainly due to their flexibility. They don't have to be completely charged every time they are used, nor do they have to be used continuously or fully drained each time.

BATEC MOBILITY only uses batteries with the safest technology available, including the BMS smart battery management system.

HANDBIKE BATEC ELECTRIC batteries don't require any maintenance.

4.12.1. INSERTING AND REMOVING THE HANDBIKE BATEC ELECTRIC BATTERY

You don't need to touch any connectors to insert or remove the battery from your HANDBIKE BATEC ELECTRIC, as this is done automatically. (See fig. 38 excerpts for Air)

To remove the battery

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Turn the key to the REMOVE BATTERY position (OPEN). Grab the handle of the battery and pull up gently until the battery comes fully out of the guide. (See fig. 39)

BATEC AIR

- Ensure the battery is turned off, with the switch in the "O" position.
- Turn the key anticlockwise by half a turn. Grip the battery, first lifting the top part and tilting it slightly, before completely removing it. (See fig. 39A)

To insert

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

Grab the battery by the handle. Set the lower right corner of the battery in the side support tab until it clicks into the guide. Slide the battery down slowly until it is resting completely in the mount. Finally, turn the key to the "OFF" position to secure the battery to the HANDBIKE BATEC ELECTRIC. (See fig. 40)

BATEC AIR

Grip the front of the battery to avoid pinching your fingers between the battery and the mount. First, place the base of the battery on the mount, ensuring the battery terminals make contact with the mount terminals. Then, attach the top part of the battery so it fits securely. Finally, turn the key clockwise by half a turn to fully secure the battery to the HANDBIKE BATEC ELECTRIC. (See fig. 40A)

INFORMATION!

On the BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT models, if the battery is in the "OFF" POSITION WITH THE KEY REMOVED, it cannot be extracted.

4.12.2. CONNECTING THE BATTERY

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

To connect the HANDBIKE BATEC ELECTRIC battery, simply turn the key to the ON position and use your HANDBIKE BATEC ELECTRIC normally. (See fig. 41)

BATEC AIR

To connect the battery to the HANDBIKE BATEC ELECTRIC, simply press the "I" power button, and use the HANDBIKE BATEC ELECTRIC as normal. (See fig. 41A)

The battery charging system is protected by the smart battery management system (BMS), which cuts power if there is a surge or if any of the battery cells are too discharged. (See fig. 42 or 42A for AIR)

CAREFUL!

Never, under any circumstances, use the discharging port on the battery to charge it or power it in any way, whether with an external power source or regenerative brake system.

4.12.3. CHARGING THE BATTERY

Recharge time for the BATEC MOBILITY battery is 3 hours. Both the charger and the battery management system cut the power when the battery is fully charged.

To charge the battery, the key must be in the OFF or OPEN position. You can charge the battery whether or not it is attached to the HANDBIKE BATEC ELECTRIC. During charging, if the battery is connected, the BATEC HANDBIKE ELECTRIC will not operate.

BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS, RAPID SPORT and SCRAMBLER SPORT

- The HANDBIKE BATEC ELECTRIC battery is charged with the 3-pin XLR port on the side of the battery. (See fig. 43)
- Open the cover protecting the charging port, connect the charger to the battery's charging port first, then plug it into the power supply and switch on the charger if it has a switch.

BATEC AIR

- The battery is charged through the coaxial port located on the top side, opposite the power switch. (See fig. 43A)
- On the top of the battery, there are LED lights that indicate the battery level when the upper button is pressed. (See fig. 43B)
- Open the cover protecting the charging port, connect the charger to the battery's charging port first, then plug it into the power supply and switch on the charger if it has a switch.

CAREFUL!

Never, under any circumstances, use chargers from other manufacturers or a charger from another HANDBIKE BATEC ELECTRIC, as the specifications for each HANDBIKE BATEC ELECTRIC series are different.

Storing the battery

Lithium ion batteries have a smart battery management system (BMS) that uses a small amount of power to monitor the battery at all times. The amount used is very small but if you store the battery fully discharged, in just a few weeks the BMS can run out of power.

If this happens, the battery can never be used again. As a result, we strongly recommend you follow these instructions carefully:

CAREFUL!

- Charge the battery fully before storing it.
- Always store the battery with the key in the OFF position.
- Charge the battery fully at least once every 3 months.
- If the battery doesn't work after being stored for a period of time, for your safety it must be disposed of.

4.12.4. BATTERY SAFETY INFORMATION

It is widely known that lithium ion batteries can be dangerous, and for good reason. Batteries used in vehicles, like those in the Batecs, have enough stored energy and potentially dangerous material to cause a significant incident if not used properly or in the case of negligence.

BATEC MOBILITY ONLY USES LITHIUM BATTERY TECHNOLOGY IT CONSIDERS STABLE AND SAFE.

All BATEC MOBILITY batteries use safety systems, including the smart battery management system (BMS), which monitors the battery at all times.

In addition to all of the built-in safety systems, to minimise the danger if a problem were to arise with elements outside the battery, such as improper functioning of the electrical installations or problems with the general power grid, we strongly recommend you follow these instructions carefully:

CAREFUL!

- Only charge the battery under supervision.
- Don't charge the battery while using your BATEC ELECTRIC HANDBIKE: Danger of electrocution and accidents.
- Don't charge the battery while asleep at night.
- Charge the battery well away from any flammable objects.
- Never use the discharging port to charge the battery and vice versa, never use the charging port to discharge it.
- Only use the BATEC MOBILITY charger provided with your HANDBIKE BATEC ELECTRIC to charge the battery.
- Never charge the battery near a heat source.
- Never charge the battery in temperatures below 0°C.
- Always store the battery fully charged.
- Don't allow the inside of the battery to be manipulated in any way, as there is a risk of explosion and/or fire.

4.12.5. ADDITIONAL INFORMATION ON THE BATTERY

Using the battery in cold conditions

The performance of lithium batteries is limited, both in terms of power and battery life, in low temperatures (below 0°C). It isn't a dangerous situation but it is bothersome, so we recommend storing the battery inside under these circumstances. Never charge the battery in temperatures below 0°C.

Using the battery in damp conditions

The casing of the BATEC MOBILITY battery provides good protection against rain and dampness for the cells and circuits inside. This protection is enough for the battery to withstand these conditions occasionally but regular, prolonged exposure to damp and rainy conditions may cause the battery to stop working properly. Thus, we recommend you avoid these circumstances and use the HANDBIKE BATEC ELECTRIC as little as possible in rainy or snowy conditions.

Battery life

The life of the battery is between 200 and 500 charging cycles, and between 1 and 3 years, depending on several factors, including temperature, frequency of use, consistency of use and other individual issues (user weight, driving style, terrain, etc.).

To get the most out of your battery we recommend you:


→ IMPORTANT!

- Drive gently at moderate speeds.
- Avoid exposure to high temperatures.
- Avoid exposure to low temperatures (below 0°C).
- Avoid leaving the battery unused for long periods of time.

The battery won't just die suddenly; it gradually begins to perform more poorly. If your battery suddenly stops working, it is normally due to a technical problem and not the end of its lifespan.

Disposing of the battery

As with all batteries, the BATEC MOBILITY batteries must be disposed of properly. Check with your national distributor or a local waste-management centre if you have any questions.

Battery indicator operation

The battery status indicator is on the display (LCD panel) of the BATEC ELECTRIC HANDBIKE. The indicator is made up of several bars that indicate the % of energy available. (See fig. 59)

→ IMPORTANT!

When the system detects that there is less than 20% battery charge available, it can limit the HANDBIKE BATEC ELECTRIC to position 1 (top speed: 6 km/h or 4 km/h, depending on the version) to save power.

4.13. WEIGHTS

Some HANDBIKE BATEC ELECTRIC models are equipped with standard weights to improve traction. The weights are located on the fork, one on each side. They each weigh 2.4 kg. To make it easier to transport the HANDBIKE BATEC ELECTRIC, the weights can be removed.

- To remove the weights:
 - 1) Hold the weight with one hand, press on the red button, tilting it so the bottom part comes away from the HANDBIKE BATEC ELECTRIC. (See fig. 44)
 - 2) Remove the weight. (See fig. 45)
- To insert the weights:
 - 1) Hold the weight with one hand. Tilt it slightly and insert the two tabs into the weight mount slot. (See fig. 46)
 - 2) Push the bottom of the weight inward until it clicks into place. (See fig. 47)

5. INSTALLATION AND ADJUSTMENTS

5.1. INSTALLATION OF THE HANDBIKE BATEC ELECTRIC ON A WHEELCHAIR

To anchor a HANDBIKE BATEC ELECTRIC to a wheelchair, you must first install a BATEC DOCKING-BAR (pat. pending). The BATEC DOCKING-BAR can be installed on most wheelchairs on the market, whether fixed or folding frame.

IMPORTANT!

An updated list of wheelchairs compatible with the HANDBIKE BATEC ELECTRIC is available in the support section of the BATEC MOBILITY website.

The HANDBIKE BATEC ELECTRIC must be installed on a wheelchair by a qualified professional authorised by BATEC MOBILITY. Contact BATEC MOBILITY to find the authorised installation professional closest to you.

5.2. USER ADJUSTMENTS

Once the HANDBIKE BATEC ELECTRIC is properly installed on the wheelchair, there are several adjustments that must be made in order to fit the HANDBIKE BATEC ELECTRIC to the user's level of mobility, type of disability and size.

These adjustments can be done easily with the tools and instructions provided with the HANDBIKE BATEC ELECTRIC and basic mechanical knowledge.

The following adjustments can be made:

- Adjusting the distance between the HANDBIKE BATEC ELECTRIC and the wheelchair.
- Adjusting the angle of the stem on the HANDBIKE BATEC ELECTRIC.
- Adjusting the height of the handlebar on the HANDBIKE BATEC ELECTRIC.
- Adjusting THE stem height
- LCD console adjustments

5.2.1. ADJUST THE DISTANCE BETWEEN THE HANDBIKE BATEC ELECTRIC AND THE CHAIR

With the HANDBIKE BATEC ELECTRIC anchored, adjust the distance between the HANDBIKE BATEC ELECTRIC and the wheelchair so that there is as little space as possible between the user's feet and the stand frame, using the telescopic anchoring bar and the screws indicated. If the user doesn't have any other preference or limitations, we recommend a distance of 8 cm. (See fig. 48 or 49A for AIR)

5.2.2. ADJUST STEM ANGLE (EXCEPT MINI 2 PLUS MODEL)

With the HANDBIKE BATEC ELECTRIC anchored, adjust the angle of the stem on the HANDBIKE BATEC ELECTRIC in order to fit the handlebar to the user's arm length, using the screws indicated. The upper and lower screws must be set on the same number. If the user doesn't have any other preference or limitations, we recommend position 3. (See fig. 49)

5.2.3. ADJUST THE HEIGHT OF THE HANDLEBAR

With the HANDBIKE BATEC ELECTRIC anchored, adjust the stem height on the HANDBIKE BATEC ELECTRIC with the Quick Release lever. (See fig. 50a or 50A for AIR)

While using the Batec, always respect the maximum height limit for the stem, not setting it above the line below the MAX. mark. It should only be taken above this line when folding the handlebar. (See fig. 50b)

CAREFUL!

Do not adjust above the MAX mark on the seat post for the BATEC MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS SCRAMBLER 2 PLUS RAPID SPORT and SCRAMBLER SPORT models.

CAREFUL!

During the height handlebar adjustment on the BATEC AIR, the handlebar cables may become trapped by components of the HANDBIKE. Ensure that this situation does not occur.

5.2.4. ADJUST THE STEM HEIGHT

MINI 2 PLUS, ELECTRIC 2 PLUS, RAPID 2 PLUS, SCRAMBLER 2 PLUS

1. With the HANDBIKE BATEC ELECTRIC unanchored, loosen the upper and lower screws connecting the fixing bar to the fork.
2. Once the fixing bar can move freely along the guide, adjust the height according to the user's chair.
3. Once the height is adjusted, re-tighten the screws so that the bar cannot move.

BATEC AIR

1. With the HANDBIKE BATEC ELECTRIC detached, loosen the upper and lower bolts securing the anchoring bar to the chassis by rotating the CLAMPING LEVERS. (See fig. 21A)

2. Once the anchoring bar can move freely along the guide, adjust the height according to the user's wheelchair, using the installation level, and secure the anchoring bar by tightening the CLAMPING LEVERS. (See fig. 21H)
3. Adjust the STOP LIMIT on the anchoring bar as shown in figure 49B by tightening the screw to 8 Nm

CAREFUL!

Before attaching the Batec Electric Air handbike, ensure that the clamping levers are securely fastened. When tightening the clamping levers, once they begin to resist, turn an additional 1/8 of a turn to fully lock the movement.

5.2.5. LCD CONSOLE ADJUSTMENTS

1. Switch on the battery to see the INFORMATION SCREEN. (See fig. 51)
2. Press the M (mode) button to move through all functions in a loop sequence from one function screen to another, such as ODO (odometer) – AVG (average speed) / MAX (maximum speed) – TRIP 1 / ODO – TRIP 2 / ODO. (See fig. 52)
 - Press R (reset) button for 2 seconds to reset TRIP 2, AVG and MAX from stored values to zero individually.
 - Data of TRIP 1, AVG and RT can be reset at the same time when one of the three data functions is being reset.
 - ODO data cannot be reset.
3. Press M and R buttons simultaneously for 2 seconds to display the SETTING SCREEN 1. (See fig. 53)
4. Press the M button to move between all functions in a loop sequence from one function screen to another, such as wheel circumference number and km/h or MPH.
 - Press R button to increment the flashing digit by 1 or convert units.
 - Press M button to confirm the digit setting and jump to next digit or next setting screen to be set.
 - The c (circumference number) values are those:
 HANDBIKE BATEC ELECTRIC MINI 2 PLUS: 990
 HANDBIKE BATEC ELECTRIC 2 PLUS / RAPID 2 PLUS / RAPID SPORT: 1445
 HANDBIKE BATEC ELECTRIC SCRAMBLER 2 PLUS / SCRAMBLER SPORT: 1621
 - Press M for 2 seconds at any setting screen to finish the setting and go to INFORMATION SCREEN
5. On SETTING SCREEN 1 (step 4) pressing M and R buttons simultaneously for 2 seconds will display the SETTING SCREEN 2. (See fig. 54)

- Battery voltage value is on the bottom and manufacturer fine adjustment battery value is on top.

6. To return to INFORMATION SCREEN and fix all settings press M button for 2 seconds, if the screen displays SETTING SCREEN 1 press M button for 2 seconds again. (See fig. 55)

IMPORTANT!

- Do not change the battery voltage value and fine adjustment battery value.
- Don't exit SETTING SCREEN 1 by switching off the battery.
- If the battery is fully charged and the screen display shows the following when switching on, you must follow steps 3, 5 and 6. (See fig. 56)

BATEC AIR

1. Turn on the battery (See fig. 56A)
2. Press the "M" button (See fig. 56B)

The display will show information on ALARM – SPEED – BATTERY LEVEL (three orange bars on the left indicate reserve; the remaining levels are green) – PARTIAL DISTANCE – TOTAL DISTANCE – ASSISTANCE LEVEL (1, 2, 3 and 4 for forward and R for reverse) – ASSISTANCE MODE (Eco, Normal, Power)

3. Screen configuration: Press both buttons  and  simultaneously for 1 second to enter the Configuration Menu (See fig. 56C).

- Use the aforementioned buttons to select the desired option from the following list:

- | | |
|------------------------------|-------------|
| 1. km/h – mph | Select |
| 2. Security pin | Disabled |
| 3. Colour | Positive |
| 4. Brightness | Maximum |
| 5. Battery info | Not enabled |
| 6. HW and SW data | |
| 7. Maximum speed | Not enabled |
| 8. Wheel circumference | 700 C |
| 9. Voltage | 36 V |
| 10. Bluetooth | Not enabled |
| 11. Reset PARTIAL DISTANCE | |
| 12. Advanced settings | Not enabled |

- Once you've selected the icon, press the "M" button to enter.
- Use the  and  buttons to select the value and confirm it by pressing the "M" button
- Use the -Select "EXIT" on the screen, then press "M" to confirm and return to the Configuration Menu.
- To return to the home screen, press both  and  simultaneously for 1 second.

6. MAINTENANCE AND REPAIRS

Your HANDBIKE BATEC ELECTRIC has been carefully designed for optimal performance, but it is important that you carry out certain preventative maintenance tasks in line with the use you give your HANDBIKE BATEC ELECTRIC. This will greatly reduce the chance of problems and increase the lifespan of your HANDBIKE BATEC ELECTRIC.

In order to minimise the risk of accidents and/or damage to your HANDBIKE BATEC ELECTRIC, it is very important to have all repairs carried out by an authorised BATEC MOBILITY distributor.

Remember that your main source of information and help is the distributor that delivered your HANDBIKE BATEC ELECTRIC, who can also respond to any questions regarding adjustments, maintenance and the warranty.

As you are a member of the "BATEC MOBILITY Family", we encourage you to stay in touch with your BATEC HANBDIKE representative in order to make the most of the product.

PREVENTATIVE MAINTENANCE PLAN SCHEDULE

TASKS	Monthly	Yearly
Check tyre pressure (2 bar). Check general tyre state (tread height, rubbing, cracking, etc.) and change if necessary.	✓	
Clean Batec.	✓	
Tighten screws. (See figs. 48, 49, 49A)	✓	
Check brakes (covers, cables, levers, discs, etc.). Change and/or adjust parts if necessary.	✓	
Check lubrication of cables and chains, grease if necessary.	✓	
Check transmission mechanism and coupling, grease if necessary.	✓	
Check spoke tension. Visit an Official Service Centre or Authorised Distributor if necessary.	✓	
Have Batec checked by an Official Service Centre or Authorised Distributor.		✓

Checking tyre pressure

Make sure the tyre pressure is as recommended for each model and that there are no irregularities on the surface.

Cleaning the HANDBIKE BATEC ELECTRIC

Clean the frame with soap and warm water. Use car soap or another non-abrasive cleaning agent. Dry it carefully. Never use a pressure washer to clean it.

Disinfecting the HANDBIKE BATEC ELECTRIC

To disinfect the HANDBIKE BATEC ELECTRIC, follow the steps below:

- Use a hydro-alcoholic sanitiser and a cloth.
- Wet the cloth with a generous amount of sanitiser.
- Rub the cloth over any surfaces that come in contact with the user: handlebars, buttons, display panel and detaching lever.
- Leave to air dry.

IMPORTANT!

Sand and salt water can damage some parts of the HANDBIKE BATEC ELECTRIC after prolonged contact. Clean it after use in these conditions.

Cleaning front and rear wheels

Make sure the front wheel and rear wheels turn freely and eliminate any dirt.

Checking brakes

Make sure that when you squeeze the brake levers, both have sufficient contact with the disc for effective braking.

Greasing cables and chains

Apply a drop of grease or light coat of Teflon spray to the cable connecting levers and to chains, after removing the old grease. After applying, remove any excess grease.

Greasing the transmission mechanism

Apply a drop of grease or light coat of Teflon spray to the transmission mechanism, after removing the old grease. After applying, remove any excess grease.

Having HANDBIKE BATEC ELECTRIC serviced by an authorised distributor

The recommended preventive maintenance and any repairs to your BATEC ELECTRIC HANDBIKE must be done by a company authorised by BATEC MOBILITY.

Replacing tyres and inner tubes

Before installing a new inner tube, check the base of the rim and inside wall of the tyre to remove any foreign bodies and, if necessary, clean well. Check the conditions of the rim tape, particularly around the valve hole. Only strong, properly placed rim tape can provide long-term protection against damage to the inner tube caused by burrs and sharp spoke heads.

Before inflating the inner tube, check the whole tube to make sure it hasn't got stuck between the tyre and the rim. Carefully pull on the valve and extract it for optimal installation.

Reusing the BATEC ELECTRIC HANDBIKE

The BATEC ELECTRIC HANDBIKE can be reused. If you would like to reuse it, please contact the distributor where you purchased it and they will take care of sending the unit to be reused to the BATEC MOBILITY facilities. The unit to be reused will be checked, given a tune-up and re-checked.

CAREFUL!

After installing the new inner tube, check tyre pressure. It is essential to keep pressure at specified levels and to make sure tyres are in proper conditions in order to ensure your personal safety and keep your HANDBIKE BATEC ELECTRIC in proper working order.

IMPORTANT!

As the necessary maintenance varies between users, check with your BATEC MOBILITY representative to establish a periodic maintenance schedule that fits your specific use.

IMPORTANT!

Always do a safety check of your HANDBIKE BATEC ELECTRIC before using. Improper maintenance may invalidate the warranty.

7. LABEL

All BATEC ELECTRIC HANDBIKES have a label on the bottom of the frame bar with the following information: (See fig. 57)

SYMBOLS

	Manufacturer
SN	Serial number
REF	Catalogue number
	This product complies with the requirements of Regulation (EU) 2017/745 on medical devices
	Date manufactured
	See user instructions
MD	With the text "Medical Device"
	Check warnings (battery)
	EN 12184:2014 – Wheelchair not intended for use as a seat in a motor vehicle
	ISO 7000 – Maximum weight allowed

When ordering spare parts or lodging claims, please provide the information on the label.

8. WARRANTY

8.1. WARRANTY PERIOD

BATEC MOBILITY guarantees all hard and structural parts for a **period of three (3) years** from the date of delivery in the case of material defects or labour for the original registered owner.

BATEC MOBILITY guarantees the BATTERY for a **period of one (1) year** from the date of delivery in the case of material defects or labour for the original registered owner. The BATTERY is considered a fungible part.

8.2. WARRANTY CONDITIONS

- Fungibles (except for the BATTERY, which has a one-year warranty) are excluded from the warranty, except in the case of premature wear on said parts caused by original manufacturing defects. These elements include tyres, inner tubes, handles, cables for brakes, gear shifter and anchoring, stand frame wheels and similar parts.
- If you have any questions regarding the warranty, BATEC MOBILITY recommends first contacting your authorised BATEC MOBILITY distributor, who will explain the solutions available for the problem.
- BATEC MOBILITY shall not be held responsible for any defects that may arise as a result of improper use or storage of the HANDBIKE BATEC ELECTRIC. If BATEC MOBILITY considers that a part has not been properly installed, used or stored, the warranty will be invalidated. Some examples are as follows:
 - The product or part hasn't been maintained according to BATEC MOBILITY recommendations, as detailed in the user instructions and maintenance instructions.
 - Non-original accessories have been used.
 - The HANDBIKE BATEC ELECTRIC or its parts have been damaged due to negligence, accident or improper use.
 - The HANDBIKE BATEC ELECTRIC or its parts have been modified in a way that is not in line with BATEC MOBILITY specifications or have undergone repairs without notifying the authorised distributor or technical service provider.
- BATEC MOBILITY reserves the right, with any claim, to require the claimant to send all parts to BATEC MOBILITY for replacement or repair at their own cost. If any of the parts are found to be defective, BATEC MOBILITY will determine whether to replace, repair or maintain the parts.

- This warranty doesn't include transport, shipping or labour costs for removing and reinstalling the parts. BATEC MOBILITY will charge the claimant for any costs incurred if the parts are found not to be defective or the damage has been caused by events not covered under this warranty.
- Defective parts will be repaired or replaced exclusively by BATEC MOBILITY or a duly qualified official BATEC MOBILITY service provider.
- BATEC MOBILITY doesn't ensure any of its products or parts will be available on the market or suitable for any particular purpose and provides no guarantee, express or implied, beyond the attached limited warranty.
- Replacement parts installed during the original warranty period shall have an additional 2-year warranty. For the BATTERY, as it is considered a fungible part, the replacement unit shall be under warranty for the remaining warranty period on the original BATTERY, up to 1 year from the date the HANDBIKE BATEC ELECTRIC was delivered.
- This warranty is non-transferable: This warranty is only valid for the original registered user, who must provide the original receipt as proof of purchase.
- BATEC MOBILITY shall not be held responsible for any products not purchased at the recommended retail price.
- Except where prohibited by law, this limited warranty constitutes the full agreement between the owner and BATEC MOBILITY regarding any sort of liability that may arise concerning any parts on the HANDBIKE BATEC ELECTRIC products.

.....

→ IMPORTANT!

- **BATEC MOBILITY** and the competent authority in the Member State of residence **should be notified of any serious incidents** related to the use of the product.
-

8.3. WARRANTY CARD

With the rest of the documentation for your HANDBIKE BATEC ELECTRIC, you will find a Warranty Card with all the details regarding your HANDBIKE BATEC ELECTRIC. You can use this card to provide all the details of your HANDBIKE BATEC ELECTRIC if a situation arises that requires the warranty.

9. REVISIONS LIST

Version	Summary / justification of changes
7.0	The RAPID SPORT and SCRAMBLER SPORT models are included.
6.0	New version BATEC ELECTRIC 2 PLUS
5.0	New Batec Air
4.0	New labelling
3.0	Updated maintenance table
2.0	Updated speeds, specifications, functions, warnings, images
1.0	To adapt to the requirements of Regulation (EU) 2017/745

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ASPHALION LTD