

The ultimate mobility tool

Batec handbikes are a comprehensive solution to your mobility problems:

- ✓ Go from reduced to **augmented mobility**; from needing help to being **independent** and **proactive**.
- ✓ Get a **2 in 1**, without giving up your manual wheelchair.
- ✓ **Overcome obstacles and drive over bumpy terrains.**

- ✓ Attach our handbikes **in just seconds** and **with no need for transfers**.
- ✓ Also available for **paraplegics**, **quadriplegics** and **hemiplegics**.



Request a quote or a demonstration without any commitment!



BATEC MOBILITY S.L.
Illa de Buda, 2
08192 Sant Quirze del Vallès
(Barcelona) - ESPAÑA
93 544 2003
info@batec-mobility.com



Batec
Electric 2 Plus
The perfect balance

Batec Electric 2 Plus

The perfect balance

The **Batec Electric 2 PLUS** offers the perfect balance with its 18-inch wheel and an exclusive tyre designed specifically for this model, **wider** and **ultra-grippy**. It's the ideal choice for both **urban environments** and **countryside trails**.



1. 18-inch wheel with an **ultra-grip tyre** and black double-wall aluminium rim.
2. **Removable QR weights** for improved traction.
3. **Anti-slip traction** control to prevent wheel spin.
4. **Cruise control** with two new dedicated buttons. Maintains a **steady speed** and disengages when braking.
5. **Parking brake**. Electronic brake with signal cut-off on both brakes.

SPECIFICATIONS:



DIAMETER:
18"



SPEED:
6km/h, 15 km/h, 20 km/h



WEIGHT:
15,9 kg (without battery)



BATTERY:
289 Wh or 576 Wh



AUTONOMY:
21,5 km* (with 280-Wh battery)
44 km* (with 576-Wh battery)

*Autonomy calculated with: 75-kg user, ultra-light wheelchair, flat terrain, no wind and a steady speed of 15 km/h. It may vary $\pm 20\%$ depending on the battery.

DIMENSIONS:

