

The ultimate mobility tool

Batec handbikes are a comprehensive solution to your mobility problems:

- ✓ Go from reduced to **augmented mobility**; from needing help to being **independent** and **proactive**.
- ✓ Get a **2 in 1**, without giving up your manual wheelchair.
- ✓ **Overcome obstacles and drive over bumpy terrains.**

- ✓ Attach our handbikes **in just seconds** and **with no need for transfers**.
- ✓ Also available for **paraplegics** and **quadriplegics**.



Request a quote or a demonstration without any commitment!



BATEC MOBILITY S.L.
Illa de Buda, 2
08192 Sant Quirze del Vallès
(Barcelona) - ESPAÑA
93 544 2003
info@batec-mobility.com



Batec Hibrid 2
Exercise for everyone

Batec Hibrid 2

Exercise for everyone

PARA QUAD
available version

The **Batec Hibrid 2** is the perfect handbike for exercising without giving up the comfort of your manual wheelchair. Together they become the **equivalent of an electric bicycle**, giving you the **extra boost** you need while you enjoy a **complete workout**.



1. **New own electronics** with 3 levels of assistance + **E mode**. Power cut-off brake sensor.
2. **Torque sensor**: Assisted starting mode on slopes without having to make a quarter turn of the pedals.
3. **New improved trigger throttle** and speed regulator.
4. **Batec Shift Lever** gear shift with 8 speeds for quadriplegic users and 9 for paraplegic users.
5. **Batec Quad Grips** to manoeuvre without taking your hands off the grips, and **Batec Magic Lever**.


 **SPEED:**
20 km/h

 **WEIGHT:**
20,6 kg (without battery)

 **AUTONOMY:**
27-40 km* (with battery 280 Wh)
55-80 km* (with battery 576 Wh)

* Autonomy calculated under the following conditions: user weighing 75 kg, ultralight wheelchair, flat terrain, no wind, and an average speed of 15 km / h. It may vary + -20% depending on the battery. The range depends on the force applied by the user.

PARA DIMENSIONS:



PARA DIMENSIONS:

