

The ultimate mobility tool

Batec handbikes are a comprehensive solution to your mobility problems:

- ✓ Go from reduced to **augmented mobility**; from needing help to being **independent** and **proactive**.
- ✓ Get a **2 in 1**, without giving up your manual wheelchair.
- ✓ **Overcome obstacles and drive over bumpy terrains.**

- ✓ Attach our handbikes **in just seconds** and **with no need for transfers**.
- ✓ Also available for **paraplegics**, **quadriplegics** and **hemiplegics**.



Request a quote or a demonstration without any commitment!



BATEC MOBILITY S.L.
Illa de Buda, 2
08192 Sant Quirze del Vallès
(Barcelona) - ESPAÑA
93 544 2003
info@batec-mobility.com





Batec
Mini 2 Plus
The very best of you

Batec Mini 2 Plus

The very best of you

Compact, practical and versatile.

The Batec **Mini 2 Plus** features the perfect combination of power, quality and ease of use in a **foldable design** that is comfortable for both transport and storage. Its **unique anchor system** means you can connect it to your wheelchair in seconds.



1. New 700-lumen front light and flashing lights on the stand frame when the **safety brake** is engaged.
2. Cruise control to maintain a **steady speed** without having to use the accelerator.
3. **Brushless motor**: 350 W nominal power, 600 W peak power.
4. **Anti-slip** traction control and **electronic braking** with signal cut-off on both brakes.
5. 12" wheel with **ultra-grip** tyre.


 **SPEED:**
20 km/h

 **WEIGHT:**
13,9 kg (without battery)

 **AUTONOMY:**
21,5 km* (With battery 280 Wh)
44 km* (With battery 576Wh)

* Autonomy calculated under the following conditions: 75 – kg user, ultra-light wheelchair, flat terrain, no wind and steady speed of 15 km/h. May vary +/- 20% depending on the battery.

FOLDED dimensions



UNFOLDED dimensions

