



Your mobility, your rules

With our **Batec handbikes**, mobility becomes **fast**, **easy** and **limitless**. In just a few seconds, you'll be on the move—no transfers, no complications. Whether you're navigating the **city** or exploring **nature**, we have a model for every adventure. **Move freely!**



The BATEC handbike made for you



More information

Number/version: MM-018-2.0 | August 2025



BATEC MOBILITY S.L.
Illa de Buda, 2 | 08192 Sant Quirze del Vallès
(Barcelona) - ESPAÑA
93 544 2003 | info@batec-mobility.com

www.batec-mobility.com



BATEC AIR, the small Batec

Batec Air,

Small, compact and ready to go


At just 11 kg and with an ultra-compact design thanks to its folding system and 10" wheel, the **Batec Air** folds in seconds. Its lightweight and practicality make it the perfect companion to take anywhere. It fits easily into its **travel bag**, while maintaining the robustness, ergonomics, and versatility you expect from Batec. **Light, agile, and always ready for action.**



TFT colour display, LED lights and buttons.




30 Lux front light.




Protective kickstand with wheels and dual rear safety lights.

DIMENSIONS

UNFOLDED



OPTIMISED FOLDING



Finished in an elegant matte black

SPECIFICATIONS



DIAMETER:
10" with aluminium rim



SPEED:
6 km/h | 10 km/h | 15 km/h



WEIGHT:
11 kg (without battery)



DISPLAY:
TFT colour



BATTERY:
216 Wh



AUTONOMY*:
23 km to 6 km/h
20 km to 10 km/h
16 km to 15 km/h



BRAKES:
• 2 mechanical disc brake callipers
• Parking brake levers
• With signal cut-off sensor

*Autonomy calculated with 75 kg user, ultralight wheelchair, flat terrain, no wind.